

4/7/99

Mike Peters  
289 Elk Creek Drive  
Bailey, CO 80421

Mike,

I have enclosed the an electronic table of the Pu and Am data, raw data for Pu and Am and copies of notebook entries Should you need anything more, please let me know

Thank you for your help

Kim Roberts  
Texas A&M University  
5007 Ave U  
Galveston, TX 77554  
(409) 741-1251

Best Available Copy

Best Available Copy

**ADMIN RECCRD**

SW-A-004058

159

May 1998.

Actinide Migration Studies at the  
Rocky Flats Environmental Technology Site.

May 1998

Second stage:

Size Distributions of Pu & Am

May 1998

## Calibration of ultrafiltration cartridges

Chemicals

System A: { 0.1  $\mu$ m  
Amicon HMP01-43  
Series #: 7174-6A3(U)

0.3  $\mu$ m Latex  
MCL

System B { 100 kDa  
Amicon S10Y100  
Series #: P7SM9/14-018

0.022  $\mu$ m  
Fluorescent Red  
Sulfate Latex

System C { 1 kDa  
Amicon S10Y1  
Series #: 12279032

Vitamin B12  
1330 Dalton

6-3-98

# Spike Ultrafiltration Experiments

## A. Labeling:

<u>SS#</u>	<u>UF System</u>	<u>Pre-filtered water</u>	<u>Permeate</u>	<u>Retentate</u>
①	0.1 $\mu$ m	PF 0.1	UF 0.1	R 0.1
②	100 kDa	PF 100	UF 100	R 100
③	1 kDa	PF 1	UF 1	R 1

## B. Experiment

Solution: tap water from Room 201 collected on June-2, 1998 at ~ 4:00 pm

Pre-filter: 0.45  $\mu$ m

Isotopes:  $^{240}\text{Pu}$  &  $^{241}\text{Am}$   
 $\sim 4 \text{ dpm}$                        $\sim 47 \text{ dpm}$

Activity:  $\sim 5 \text{ dpm}$  for each isotope;

6-3-98.

<u>SS#</u>	<u>System</u>	<u>Vol. of PFW (L)</u>	<u>Vol. of Coal</u>	<u>Starting Time/End</u>	<u>Pressure (PSI)</u>
(1)	0.1 $\mu$ m H.M.Pol-43	~20	2.25	3:00 pm / 5:10 pm	straight valve
(2)	100 kDa (SiO <sub>2</sub> 100)	~20	~2.10	2:40 pm / 3:40 pm	straight valve
(3)	1 kDa (SiO <sub>2</sub> 1)	~20	2.25	2:50 pm / 4:10 pm	~30 psi

~~102~~ Walnut Cr. ABV  
Flume Pond @ 6505

8/26/98.  
CSM

A: SPM

<u>SS#</u>	<u>Filter #</u>	<u>Water Vol.</u>	C-G
1	E-011	$300 + (-70ml) = 230ml$	$\div 10 = 240ml$
2	E-012	$105ml$	$\div 10 = 105ml$ only
3	E-013		

B: POC (filter before 1992)

<u>SS#</u>	<u>Filter #</u>	<u>Vol. (ml)</u>
1	27	$100 + \rightarrow 100ml$
2	28	52ml
3	26	BIK

C: DOC

<u>SS#</u>	<u>Bottle #</u>	<u>Acidify</u>
1	B 275	✓
2	B 062	✓
3	C 068	✓
4		

pH. = 9.48

(2.5 by pH paper)

DO<sub>2</sub> = 7.72 mg/L

Water Temp = 25°C

Air Temp = 31.5°C

Specific Conductivity = 0.556 mS/cm

Co. 45  $\mu$ m @ GS05

8-26-98.

→ 0.1  $\mu$ m System - & 100 kDa System  
20 l start @ 4:30 pm stop @ 3:00 pm (8-27/98) @ 3.65 l

→ 100 kDa Cartridge - II { start @ 8:03 pm  
end @ 9:45 pm }

Combined coll  
&  
processed  
further to  
3.4 l

→ 100 kDa Cartridge - I. (20 l  $\rightarrow$  3.2 l Ret.)

SS#	Bottle #	Fraction	Acid
1	A 527	UFW-100	✓ 3d
2	B 291	UFW-100	✓
3	B 233	COC100-1 @ 3.2 l	✓ 2d
4	B 032	COC100-2 @ 3.2 l	✓ 2d

Final concentration Factor for 100 kDa System -  
A. 100 kDa: ~40 l  $\rightarrow$  3.4 liter of Coll.

B. 0.1  $\mu$ m System ~20 l  $\rightarrow$  3.65 l

SS#	DOC Bottle #	Fraction	Acid
1	A 161	UFW 0.1 20 ml	✓ 3d
2	D 001	UF 0.1 "	✓ "
3	B 016	COC 0.1 15 ml	✓ 2d
4	A 599	COC 0.1	✓ ✓ "
	A 599		

8-26-78  
G503 site.

① Particle & filtration.

$$\frac{2.6 \text{ Gal} + 1.6 \text{ Gal}}{2} = \overset{2.10}{\sim 2.1} \text{ Gal} \approx 7.8 \text{ l}$$

Volume for  $\geq 20 \mu\text{m}$  filters. (x3 filters)  
⊗  $\geq 0.45 \mu\text{m}$  particle loss a little bit in #2 & 3 filters.

② Total water for ppt:  $\sim 3.5 \text{ l}$

③ water ( $< 0.45 \mu\text{m}$ ) for ppt  
 $\leq 20$  liters

④

8-27-98.

A: SPM

<u>SS #</u>	<u>Filter #</u>	<u>Vol. (ml)</u>	<u>Mark</u>
1	E-013	100 +	✓
2	E-014	100 ml	✓
3	E-015	100 ml	✓
4	E-016	100 ml	✓

B: POC

<u>SS #</u>	<u>Filter #</u>	<u>Vol. (ml)</u>	<u>Acid</u>	<u>Mark</u>
1	<1992-29	50 ml ✓	✓	98827-POC-01
2	-30	50 ml ✓	-	98827-POC-02
3	-31	50 ml ✓	✓	POC-03
4	-32	58 ml (58 ml) ✓	-	POC-04
5/6	-33/34	Blk		

C: DOC

<u>SS #</u>	<u>Bottle #</u>	<u>Vol (ml)</u>	<u>Acid/FA</u>
1	B153	20 ml	✓
2	B059	"	"
3	A185	"	"
4	B054	"	"

100 kDa

98827A.

1998-8-27

Prefiltration: mixed up & prefilter it (5um → 0.45um)  
into 8 20L containers

- i) #1 & #8 for 98827A 100 kDa.
- ii) #2 & #7 for 98827B 100 kDa;
- iii) #3 & #6 for <0.45um total; TD (total Diss.)
- iv) #4 for 98827A 0.1um system;
- v) #5 for 98827B 0.1um system;

100 kDa 98827A

start @ 5:10 pm → 6:45 pm

{ i) 20L → 1.80L  
 ii) 20L → ≥1.8L (≥1.8L)

Sample for DOC }

start @ 7:07 pm → 8:34 pm

100L → 4.2L { 98827B . i) start @ { 9:50am → 11:20am } ✓ Sample for DOC  
 { 20L → 2.35L  
 ii) start @ { start @ 12:50pm → 2:15pm  
 { 20L → 1.8L

0.1um System 98827A. #4 PFW

i) start @ 8:00 pm → 7:00 am (8-28-98)  
 20L → 3.2L (3.2L)

ii) 0.1um System 98827B #5 PFW

start @ 1:35 pm → 4:05 AM (8-29-98)  
 20L → ≥2.6L

DOC from UFZ

8-27-98

98827A

<u>SS#</u>	<u>Bottle #</u>	<u>Fraction</u>	<u>Acid.</u>	<u>Mark.</u>
1	C190	100 U <sub>1</sub> - <del>A</del>	✓	A don't
2	C191	100 U <sub>1</sub> -II	✓	-
3	B127	100 R <sub>1</sub> -I	✓	15ml
4	C068	100 R <sub>1</sub> -II	✓	"
5	C006	100 U <sub>2</sub>	✓	20 ml
6	B191	100 R <sub>2</sub>	✓	15 ml
2	B124	0.14 U <sub>1</sub>	(98827A) ✓ <sup>2d</sup>	20 ml
3	A201	0.14 R <sub>1</sub>	(98827A) ✓ <sup>2d</sup>	15 ml

~~9~~  
~~10~~


98827B

9	D020	B100 U <sub>1</sub>	✓	20 ml
10	B119	B100 R <sub>1</sub>	✓	15 ml
11	A105	B0.1 U <sub>1</sub> (UFW)	✓	20 ml
12	B286	B0.1 R <sub>1</sub> (coll.)	✓	15 ml


8-27-90

total water filtered

98827A



98827B



Solid Sample Digestion (SOP 773)

- 1 Combust sample if applicable
- 2 Add tracers
- 3 In teflon beaker add 20 ml conc HNO<sub>3</sub>, 20 ml conc HCL and 20 ml conc HF
- 4 Ensure solids go in to solution (IF not repeat step 3 ) Evaporate to dryness
- 5 Add 10 ml conc HNO<sub>3</sub> and 5 ml 30% H<sub>2</sub>O<sub>2</sub> Let stand at room temp for a few minutes then evaporate to dryness
- 6 Add 10 ml conc HNO<sub>3</sub>, 100ml DI and 2g H<sub>3</sub>BO<sub>3</sub> Heat to ensure dissolution
- 7 Transfer to 250 ml cent Tube using DI
- 8 Add 1 ml (40mg/ml) Fe carrier
- 9 Ppt Fe(OH)<sub>3</sub> by adding conc NH<sub>4</sub>OH (~50ml)
- 10 Centrifuge at 3500 rpm for 30 minutes
- 11 Dissolve ppt With 3x the vol of ppt with conc HCL
- 12 Add 75 ml 9N HCL
- 13 Add 1 ml saturated NaNO<sub>2</sub> mix thoroughly

	C98826 5 (RF4)	C98826 0.45 (RF6)	C98827 5 (RF10)	C98827 0.45 (RF18)	Blankspk4 (RF11)	Blank (RF12)
	9/6/18					→
	9/7/18					→
	✓	✓	✓	✓	✓	✓
	✗	✓	✗	✓	✓	✓
		9/8/18		9/8/18		→
		✓		✓	✓	✓
		✓		✓	✓	✓
		✓		✓	✓	✓
		✓		✓	✓	✓
		7ml 9/9/18		5ml 9/9/18	2ml	2ml
		✓		✓	✓	✓
		✓		✓	✓	✓

✗ repeated digestion w/ conc HCl & HNO<sub>3</sub> moved to bath of other solids to start digestion from beginning.

**Pu column (SOP 778)**

- 1 Fill disposable plastic column with 7 cm AG1x8 resin (by resin/DI slurry)
- 2 Place funnel on top of column with Whatman filter paper
- 3 Wet filter paper with 9N HCl
- 4 Condition resin with 50 ml 9 N HCl
- 5 Load sample through filter, rinse 2x with 20 ml 9N HCl (= Am fraction)
- 6 Rinse 2x with 20 ml 9N HCl, discard
- 7 Elute Pu with 20ml 9N HCl + 1.5 ml HI
- 8 Add 1 ml conc HNO<sub>3</sub>, evaporate to dryness

	RF 4	RF 6	RF 16	RF 18	RF 11	RF 12
	<del>9/9/8</del>	9/9/8	<del>9/9/8</del>	9/9/8		
		✓		✓	✓	✓
		✓		✓	✓	✓
		✓		✓	✓	✓
		✓		✓	✓	✓
		✓		✓	✓	✓
		✓		✓	✓	✓
		✓		✓	✓	✓

**Pu Microprecipitation (SOP 778)**

- 1 Add 1 ml conc HCl, mix well to resolubilize
- 2 Add 14 ml DI, mix well
- 3 Add 1.0 ml lanthanum carrier + 0.5 ml H<sub>2</sub>O<sub>2</sub> mix well
- 4 Add 5 ml 3N HF
- 5 Let sample sit for min 15-20 minutes
- 6 Set up filtration apparatus place 25 mm filter membrane on support screen and lock
- 7 Apply vacuum, rinse filter with 1-2 ml methanol, rinse filter with DI
- 8 Transfer sample, rinsing beaker once with 5 ml DI
- 9 Rinse filter with 10-15 ml DI
- 10 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 11 Mount filter on stainless steel planchet with double sided adhesive tape

	✓		✓	✓	✓	
	✓		✓	✓	✓	
	✓		✓	✓	✓	
	✓		✓	✓	✓	
	✓		✓	✓	✓	
	✓		✓	✓	✓	
	✓		✓	✓	✓	
	✓		✓	✓	✓	
	✓		✓	✓	✓	

**Am Methanolic Anion Exchange Column (SOP 780)**

- 1 Mix anion exchange resin with twice the volume of 1N HNO<sub>3</sub>/93% methanol solution (for 160 ml 10ml conc HNO<sub>3</sub>, 150 ml methanol) overnight
- 2 Add 5-10 ml conc HNO<sub>3</sub>,
- 3 Add 100-125 ml DI,
- 4 Add 10 ml Fe carrier
- 5 Transfer to cent tube and ppt with conc NH<sub>4</sub>OH Cent @3500 rpm for 15 min/
- 6 Add conc HNO<sub>3</sub> and evaporate, Repeat
- 7 Add conc HNO<sub>3</sub> until dissolved Add methanol (15 ml of methanol for each 1 ml conc HNO<sub>3</sub>)
- 8 Pour resin slurry into disposable column and let resin settle to 7cm, place a layer of silica or glass beads on top of resin
- 9 Place funnel on top of column with Whatman filter paper, wet filter with 1N HNO<sub>3</sub> /93% methanol solution
- 10 Condition column with 40 ml of 1N HNO<sub>3</sub> /93% methanol solution
- 8 Load the sample onto column through filter, rinse with 25 ml 1N HNO<sub>3</sub> /93% methanol solution
- 9 Remove and discard filter and rinse twice with 25 ml 1N HNO<sub>3</sub> / 93% methanol solution Discard into methanolic waste receptacle
- 10 Strip Am by passing three 20 ml volumes of 8 N HNO<sub>3</sub> through the column allowing each rinse to pass completely before adding next rinse Collect eluate in beaker for Teva column
- 11 Evaporate to dryness

	RF4	RF9	RF10	RF18	RF11	RF12
		9/10/98		9/10/98		→
		✓		✓	✓	✓
		✓		✓	✓	✓
		✓		✓	✓	✓
		✓		✓	✓	✓
		9/11/98		9/11/98		→
		9/12/98		9/12/98		→
		✓		✓	✓	✓
		✓		✓	✓	✓
		✓		✓	✓	✓
		"		"	"	"
		<del>9/12/98</del>		<del>9/12/98</del>	9/13/98	9/14/98

**Am Teva Resin (SOP 780)**

- 1 Redissolve sample from above in 10 ml 2M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100ml 15 2 g NH<sub>4</sub>SCN & 0 35 ml 98% formic acid in 100ml DI) Allow sample to sit for 1 hour to ensure dissolution
- 2 Condition a TEVA resin 2 ml column with 5 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid solution
- 3 Transfer sample into TEVA column in two portions using disposable polyethylene transfer pipet Rinse the sample container with 1 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid and transfer to column Repeat rinse and add to column
- 4 Rinse the TEVA column with two 5 ml volumes of 1M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100 ml dissolve 7 6 g NH<sub>4</sub>SCN and 0 35 ml 98% formic acid in 100 ml of DI) Allow first wash to pass completely before adding the second wash (This washes lanthanides from column)
- 5 Strip Am from column with 15 ml of 2N HCl in three 5ml portions allowing each 5 ml to pass completely
- 6 To decompose thiocyanate, add 2 5 ml conc HNO<sub>3</sub> and 7 5 ml conc HCL to the Am solution Swirl gently Evaporate until ~1 drop solution is remaining
- 7 Add 5 ml conc HNO<sub>3</sub> Evaporate until volume is ~1 drop

	<del>RF7</del>	RF9	<del>RF10</del>	RF18	RF11	RF12
	<del>✓</del>	9/13/98	<del>✓</del>	9/13/98		→
	<del>✓</del>	9/14/98	<del>✓</del>	9/14/98		→
	<del>✓</del>	✓	<del>✓</del>	✓	✓	✓
	<del>✓</del>	"	<del>✓</del>	"	"	"
	<del>✓</del>		<del>✓</del>			
	<del>✓</del>	"	<del>✓</del>	"	"	"
	<del>✓</del>	✓	<del>✓</del>	✓	✓	✓

**Am Micro-precipitation (SOP 780)**

- 1 Add 1 ml conc HCl to sample Heat for 5 minutes Add 15 ml DI
- 2 Add 0 5 ml lanthanum carrier mix well Add 5 ml HF Mix well
- 3 Allow sample to stand for 15-20 minutes minimum
- 4 Place 25 mm filter membrane in a filter funnel assembly and turn on vacuum Rinse with 1-2 ml alcohol
- 5 Load sample into filter Rinse sample beaker once with 5 ml DI and add to funnel
- 6 After sample has passed through filter Rinse filter with 10-15 ml DI
- 7 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 8 Mount filter on stainless steel planchet with double sided adhesive tape

	RF9	RF18	RF11	RF12
	✓	✓	✓	✓
	✓	✓	✓	✓
	✓	✓	✓	✓
	✓	✓	✓	✓
	✓	✓	✓	✓
	✓	✓	✓	✓
	✓	✓	✓	✓

Completed - 9/15/98

**Solid Sample Digestion (SOP 773)**

- 1 Combust sample if applicable *leached in 8N HNO<sub>3</sub>*
- 2 Add tracers
- 3 In teflon beaker add 20 ml conc HNO<sub>3</sub>, 20 ml conc HCL and 20 ml conc HF
- 4 Ensure solids go in to solution (IF not repeat step 3 ) Evaporate to dryness
- 5 Add 10 ml conc HNO<sub>3</sub> and 5 ml 30% H<sub>2</sub>O<sub>2</sub> Let stand at room temp for a few minutes then evaporate to dryness
- 6 Add 10 ml conc HNO<sub>3</sub>, 100ml DI and 2g H<sub>3</sub>BO<sub>3</sub> Heat to ensure dissolution
- 7 Transfer to 250 ml cent Tube using DI
- 8 Add 1 ml (40mg/ml) Fe carrier
- 9 Ppt Fe(OH)<sub>3</sub> by adding conc NH<sub>4</sub>OH (~50ml)
- 10 Centrifuge at 3500 rpm for 30 minutes
- 11 Dissolve ppt With 3x the vol of ppt with conc HCL
- 12 Add 75 ml 9N HCl
- 13 Add 1 ml saturated NaNO<sub>2</sub> mix thoroughly

	F98826 7045 (RF5)	F98826 720 (RF3)	F98827A 7045 (RF17)	F98827A 720 (RF15)	F98827B 720 (RF26)
9/6/98					↓
9/7/98					↓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
9/8/98					→
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	1.5ml	5ml	2ml	3ml	3ml
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓

**Pu column (SOP 778)**

- 1 Fill disposable plastic column with 7 cm AG1x8 resin (by resin/DI slurry)
- 2 Place funnel on top of column with Whatman filter paper
- 3 Wet filter paper with 9N HCl
- 4 Condition resin with 50 ml 9 N HCl
- 5 Load sample through filter, rinse 2x with 20 ml 9N HCl (= Am fraction)
- 6 Rinse 2x with 20 ml 9N HCl, discard
- 7 Elute Pu with 20ml 9N HCl + 1.5 ml HI
- 8 Add 1 ml conc HNO<sub>3</sub>, evaporate to dryness

	RF5	RF3	RF17	RF15	RF20	
9/9/98					→	
✓	✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	✓	

**Pu Microprecipitation (SOP 778)**

- 1 Add 1 ml conc HCl, mix well to resolubilize
- 2 Add 14 ml DI, mix well
- 3 Add 1.0 ml lanthanum carrier + 0.5 ml H<sub>2</sub>O<sub>2</sub> mix well
- 4 Add 5 ml 3N HF
- 5 Let sample sit for min 15-20 minutes
- 6 Set up filtration apparatus place 25 mm filter membrane on support screen and lock
- 7 Apply vacuum, rinse filter with 1-2 ml methanol, rinse filter with DI
- 8 Transfer sample, rinsing beaker once with 5 ml DI
- 9 Rinse filter with 10-15 ml DI
- 10 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 11 Mount filter on stainless steel planchet with double sided adhesive tape

✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	
	✓	✓	✓	✓	
	✓	✓	✓	✓	
	✓	✓	✓	✓	

**Am Methanolic Anion Exchange Column (SOP 780)**

- 1 Mix anion exchange resin with twice the volume of 1N HNO<sub>3</sub>/93% methanol solution (for 160 ml 10ml conc HNO<sub>3</sub>, 150 ml methanol) overnight
- 2 Add 5-10 ml conc HNO<sub>3</sub>,
- 3 Add 100-125 ml DI,
- 4 Add 10 ml Fe carrier
- 5 Transfer to cent tube and ppt with conc NH<sub>4</sub>OH Cent @3500 rpm for 15 min/
- 6 Add conc HNO<sub>3</sub> and evaporate, Repeat
- 7 Add conc HNO<sub>3</sub> until dissolved Add methanol (15 ml of methanol for each 1 ml conc HNO<sub>3</sub>)
- 8 Pour resin slurry into disposable column and let resin settle to 7cm, place a layer of silica or glass beads on top of resin
- 9 Place funnel on top of column with Whatman filter paper, wet filter with 1N HNO<sub>3</sub> /93% methanol solution
- 10 Condition column with 40 ml of 1N HNO<sub>3</sub> /93% methanol solution
- 8 Load the sample onto column through filter, rinse with 25 ml 1N HNO<sub>3</sub> /93% methanol solution
- 9 Remove and discard filter and rinse twice with 25 ml 1N HNO<sub>3</sub> / 93% methanol solution Discard into methanolic waste receptacle
- 10 Strip Am by passing three 20 ml volumes of 8 N HNO<sub>3</sub> through the column allowing each rinse to pass completely before adding next rinse Collect eluate in beaker for Teva column
- 11 Evaporate to dryness

	RF5	RF3	RF17	RF15	RF26
9/10/98	→				
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
9/11/98	→				
9/12/98	→				
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
9/13/98	9/13/98	9/12/98	9/12/98	9/13/98	9/13/98

**Am Teva Resin (SOP 780)**

- 1 Redissolve sample from above in 10 ml 2M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100ml 15 2 g NH<sub>4</sub>SCN & 0 35 ml 98% formic acid in 100ml DI) Allow sample to sit for 1 hour to ensure dissolution
- 2 Condition a TEVA resin 2 ml column with 5 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid solution
- 3 Transfer sample into TEVA column in two portions using disposable polyethylene transfer pipet Rinse the sample container with 1 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid and transfer to column Repeat rinse and add to column
- 4 Rinse the TEVA column with two 5 ml volumes of 1M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100 ml dissolve 7 6 g NH<sub>4</sub>SCN and 0 35 ml 98% formic acid in 100 ml of DI) Allow first wash to pass completely before adding the second wash (This washes lanthanides from column)
- 5 Strip Am from column with 15 ml of 2N HCl in three 5ml portions allowing each 5 ml to pass completely
- 6 To decompose thiocyanate, add 2 5 ml conc HNO<sub>3</sub> and 7 5 ml conc HCL to the Am solution Swirl gently Evaporate until ~1 drop solution is remaining
- 7 Add 5 ml conc HNO<sub>3</sub> Evaporate until volume is ~1 drop

	RF5	RF3	RF7	RF15	RF20
	9/13/98				→
	✓	✓	✓	✓	✓
	9/14/98				→
	✓	✓	✓	✓	✓
	"	"	"	"	"
	"	"	"	"	"
	✓	✓	✓	✓	✓

**Am Micro-precipitation (SOP 780)**

- 1 Add 1 ml conc HCl to sample Heat for 5 minutes Add 15 ml DI
- 2 Add 0 5 ml lanthanum carrier mix well Add 5 ml HF Mix well
- 3 Allow sample to stand for 15-20 minutes minimum
- 4 Place 25 mm filter membrane in a filter funnel assembly and turn on vacuum Rinse with 1-2 ml alcohol
- 5 Load sample into filter Rinse sample beaker once with 5 ml DI and add to funnel
- 6 After sample has passed through filter Rinse filter with 10-15 ml DI
- 7 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 8 Mount filter on stainless steel planchet with double sided adhesive tape

	RF5	RF3	RF7	RF15	RF20
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓

Completed 9/15/98



**Pu column (SOP 778)**

- 1 Fill disposable plastic column with 7 cm AG1x8 resin (by resin/DI slurry)
- 2 Place funnel on top of column with Whatman filter paper
- 3 Wet filter paper with 9N HCl
- 4 Condition resin with 50 ml 9 N HCl
- 5 Load sample through filter, rinse 2x with 20 ml 9N HCl (= Am fraction)
- 6 Rinse 2x with 20 ml 9N HCl, discard
- 7 Elute Pu with 20ml 9N HCl + 1.5 ml HI
- 8 Add 1 ml conc HNO<sub>3</sub>, evaporate to dryness

	RF8	RF10	RF20	RF21	RF22
1	✓	✓	✓	✓	✓
2	✓	✓	✓	✓	✓
3	✓	✓	✓	✓	✓
4	✓	✓	✓	✓	✓
5	✓	✓	✓	✓	✓
6	✓	✓	✓	✓	✓
7	✓	✓	✓	✓	✓
8	✓	✓	✓	✓	✓

**Pu Microprecipitation (SOP 778)**

- 1 Add 1 ml conc HCl, mix well to resolubilize
- 2 Add 14 ml DI, mix well
- 3 Add 1.0 ml lanthanum carrier + 0.5 ml H<sub>2</sub>O<sub>2</sub>, mix well
- 4 Add 5 ml 3N HF
- 5 Let sample sit for min 15-20 minutes
- 6 Set up filtration apparatus place 25 mm filter membrane on support screen and lock
- 7 Apply vacuum, rinse filter with 1-2 ml methanol, rinse filter with DI
- 8 Transfer sample, rinsing beaker once with 5 ml DI
- 9 Rinse filter with 10-15 ml DI
- 10 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 11 Mount filter on stainless steel planchet with double sided adhesive tape

	RF8	RF10	RF20	RF21	RF22
1	✓	✓	✓	✓	✓
2	✓	✓	✓	✓	✓
3	✓	✓	✓	✓	✓
4	✓	✓	✓	✓	✓
5	✓	✓	✓	✓	✓
6	✓	✓	✓	✓	✓
7	✓	✓	✓	✓	✓
8	✓	✓	✓	✓	✓
9	✓	✓	✓	✓	✓
10	✓	✓	✓	✓	✓
11	✓	✓	✓	✓	✓

**Am Methanolic Anion Exchange Column (SOP 780)**

- 1 Mix anion exchange resin with twice the volume of 1N HNO<sub>3</sub>/93% methanol solution (for 160 ml 10ml conc HNO<sub>3</sub>, 150 ml methanol) overnight
- 2 Add 5-10 ml conc HNO<sub>3</sub>,
- 3 Add 100-125 ml DI,
- 4 Add 10 ml Fe carrier
- 5 Transfer to cent tube and ppt with conc NH<sub>4</sub>OH Cent @3500 rpm for 15 min/
- 6 Add conc HNO<sub>3</sub> and evaporate, Repeat
- 7 Add conc HNO<sub>3</sub> until dissolved Add methanol (15 ml of methanol for each 1 ml conc HNO<sub>3</sub>)
- 8 Pour resin slurry into disposable column and let resin settle to 7cm, place a layer of silica or glass beads on top of resin
- 9 Place funnel on top of column with Whatman filter paper, wet filter with 1N HNO<sub>3</sub> /93% methanol solution
- 10 Condition column with 40 ml of 1N HNO<sub>3</sub> /93% methanol solution
- 8 Load the sample onto column through filter, rinse with 25 ml 1N HNO<sub>3</sub> /93% methanol solution
- 9 Remove and discard filter and rinse twice with 25 ml 1N HNO<sub>3</sub> / 93% methanol solution Discard into methanolic waste receptacle
- 10 Strip Am by passing three 20 ml volumes of 8 N HNO<sub>3</sub> through the column allowing each rinse to pass completely before adding next rinse Collect eluate in beaker for Teva column
- 11 Evaporate to dryness

	RF8	RF10	RF20	RF21	RF22
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	<del>7/17/98</del> →				
	✓	✓	✓	✓	✓
	<del>9/21/98</del> →				
	<del>9/23/98</del> →				
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	<del>9/24/98</del> →				

**Am Teva Resin (SOP 780)**

- 1 Redissolve sample from above in 10 ml 2M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100ml 15 2 g NH<sub>4</sub>SCN & 0 35 ml 98% formic acid in 100ml DI) Allow sample to sit for 1 hour to ensure dissolution
- 2 Condition a TEVA resin 2 ml column with 5 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid solution
- 3 Transfer sample into TEVA column in two portions using disposable polyethylene transfer pipet Rinse the sample container with 1 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid and transfer to column Repeat rinse and add to column
- 4 Rinse the TEVA column with two 5 ml volumes of 1M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100 ml dissolve 7 6 g NH<sub>4</sub>SCN and 0 35 ml 98% formic acid in 100 ml of DI) Allow first wash to pass completely before adding the second wash (This washes lanthanides from column)
- 5 Strip Am from column with 15 ml of 2N HCl in three 5ml portions allowing each 5 ml to pass completely
- 6 To decompose thiocyanate, add 2 5 ml conc HNO<sub>3</sub> and 7 5 ml conc HCL to the Am solution Swirl gently Evaporate until ~1 drop solution is remaining
- 7 Add 5 ml conc HNO<sub>3</sub> Evaporate until volume is ~1 drop

RF8	RF10	RF20	RF21	RF22
9/24/98				✓
✓	✓	✓	✓	
✓	✓	✓	✓	
✓	✓	✓	✓	
9/25/98			→	

**Am Micro-precipitation (SOP 780)**

- 1 Add 1 ml conc HCl to sample Heat for 5 minutes Add 15 ml DI
- 2 Add 0 5 ml lanthanum carrier mix well Add 5 ml HF Mix well
- 3 Allow sample to stand for 15-20 minutes minimum
- 4 Place 25 mm filter membrane in a filter funnel assembly and turn on vacuum Rinse with 1-2 ml alcohol
- 5 Load sample into filter Rinse sample beaker once with 5 ml DI and add to funnel
- 6 After sample has passed through filter Rinse filter with 10-15 ml DI
- 7 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 8 Mount filter on stainless steel planchet with double sided adhesive tape

RF8	RF10	RF20	RF21	RF22
9/25/98				✓
✓	✓	✓	✓	✓
✓	✓	✓	✓	✓
✓	✓	✓	✓	✓
✓	✓	✓	✓	✓
✓	✓	✓	✓	✓
✓	✓	✓	✓	✓

⊗ filtered v. quickly

## Rocky Flats Am and Pu Sample Procedures

### Coprecipitation (SOP 778)

- 1 To 10L of sample add of conc HNO<sub>3</sub> until pH < 2 Let sample sit for 16 hours after acidification before processing
- 2 Add 50 ml conc HCl
- 3 Add tracers (<sup>243</sup>Am and <sup>242</sup>Pu)
- 4 Add 5 ml ferric chloride carrier (40 mg Fe<sup>3+</sup>/ml)
- 5 Mix sample and measure pH with pH paper, if pH > 1 add 12N HCl until pH < 1 Cover and stir for 30 min
- 6 Check pH again If < 1, gently add conc NH<sub>4</sub>OH, add until turbidity persists than add additional 50 ml
- 7 Stir for 30 minutes then let settle overnight
- 8 Siphon ~9L of solution, discard
- 9 Transfer remaining ~1L into 250ml centrifuge tubes, rinsing bucket and transfer beaker (if used) with minimum DI
- 10 Centrifuge for 15 min at 3500 rpm
- 11 Decant and discard supernate
- 12 To dissolve precipitate add 3 times ppt volume of conc HCl and mix by vortexing (Shipped in 3N HCl → evaporated)
- 13 Add 75 ml of 9 N HCl
- 14 Add 2 ml of saturated NaNO<sub>2</sub> to samples, mix well and set aside for 15 minutes

	RF2	RF7	RF9	RF14	RF19	RF27	RF31	RF32	
	9/2/98 →								
	✓	✓	✓	✓	✓	✓	✓	✓	
	✓	✓	✓	✓	✓	✓	✓	✓	

**Pu column (SOP 778)**

- 1 Fill disposable plastic column with 7 cm AG1x8 resin (by resin/DI slurry)
- 2 Place funnel on top of column with Whatman filter paper
- 3 Wet filter paper with 9N HCl
- 4 Condition resin with 50 ml 9 N HCl
- 5 Load sample through filter, rinse 2x with 20 ml 9N HCl (= Am fraction)
- 6 Rinse 2x with 20 ml 9N HCl, discard
- 7 Elute Pu with 20ml 9N HCl + 1.5 ml HI
- 8 Add 1 ml conc HNO<sub>3</sub>, evaporate to dryness

RF2	RF7	RF9	RF14	RF19	RF27	RF31	RF32
9/11/98							
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓

**Pu Microprecipitation (SOP 778)**

- 1 Add 1 ml conc HCl, mix well to resolubilize
- 2 Add 14 ml DI, mix well
- 3 Add 1.0 ml lanthanum carrier + 0.5 ml H<sub>2</sub>O<sub>2</sub>, mix well
- 4 Add 5 ml 3N HF
- 5 Let sample sit for min 15-20 minutes
- 6 Set up filtration apparatus place 25 mm filter membrane on support screen and lock
- 7 Apply vacuum, rinse filter with 1-2 ml methanol, rinse filter with DI
- 8 Transfer sample, rinsing beaker once with 5 ml DI
- 9 Rinse filter with 10-15 ml DI
- 10 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 11 Mount filter on stainless steel planchet with double sided adhesive tape

RF2	RF7	RF9	RF14	RF19	RF27	RF31	RF32
✓	9/12/98	✓	✓	✓	✓	9/14/98	✓
✓	✓	✓	✓	✓	✓	✓	✓
9/15/98							
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓

Am Methanolic Anion Exchange Column (SOP 780)

- 1 Mix anion exchange resin with twice the volume of 1N HNO<sub>3</sub>/93% methanol solution (for 160 ml 10ml conc HNO<sub>3</sub>, 150 ml methanol) overnight
- 2 Add 5-10 ml conc HNO<sub>3</sub>,
- 3 Add 100-125 ml DI,
- 4 Add 10 ml Fe carrier
- 5 Transfer to cent tube and ppt with conc NH<sub>4</sub>OH Cent @3500 rpm for 15 min/
- 6 Add conc HNO<sub>3</sub> and evaporate, Repeat
- 7 Add conc HNO<sub>3</sub> until dissolved Add methanol (15 ml of methanol for each 1 ml conc HNO<sub>3</sub>)
- 8 Pour resin slurry into disposable column and let resin settle to 7cm, place a layer of silica or glass beads on top of resin
- 9 Place funnel on top of column with Whatman filter paper, wet filter with 1N HNO<sub>3</sub> /93% methanol solution
- 10 Condition column with 40 ml of 1N HNO<sub>3</sub> /93% methanol solution
- 8 Load the sample onto column through filter, rinse with 25 ml 1N HNO<sub>3</sub> /93% methanol solution
- 9 Remove and discard filter and rinse twice with 25 ml 1N HNO<sub>3</sub> / 93% methanol solution Discard into methanolic waste receptacle
- 10 Strip Am by passing three 20 ml volumes of 8 N HNO<sub>3</sub> through the column allowing each rinse to pass completely before adding next rinse Collect eluate in beaker for Teva column
- 11 Evaporate to dryness

	RF2	RF7	RF9	RF14	RF19	RF27	RF31	RF32	
11/2/98	→						9/2/98		
✓	✓	✓	✓	✓	✓		✓	✓	
✓	✓	✓	✓	✓	✓		✓	✓	
9/14/98	→							→	
✓	✓	✓	✓	✓	✓		✓	✓	
7/16/98	→								→
✓	✓	✓	✓	✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	✓	✓	✓	✓	
✓	✓	✓	✓	✓	✓	✓	✓	✓	

**Am Teva Resin (SOP 780)**

- 1 Redissolve sample from above in 10 ml 2M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100ml 15 2 g NH<sub>4</sub>SCN & 0 35 ml 98% formic acid in 100ml DI) Allow sample to sit for 1 hour to ensure dissolution
- 2 Condition a TEVA resin 2 ml column with 5 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid solution
- 3 Transfer sample into TEVA column in two portions using disposable polyethylene transfer pipet Rinse the sample container with 1 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid and transfer to column Repeat rinse and add to column
- 4 Rinse the TEVA column with two 5 ml volumes of 1M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100 ml dissolve 7 6 g NH<sub>4</sub>SCN and 0 35 ml 98% formic acid in 100 ml of DI) Allow first wash to pass completely before adding the second wash (This washes lanthanides from column)
- 5 Strip Am from column with 15 ml of 2N HCl in three 5ml portions allowing each 5 ml to pass completely
- 6 To decompose thiocyanate, add 2 5 ml conc HNO<sub>3</sub> and 7 5 ml conc HCL to the Am solution Swirl gently Evaporate until ~1 drop solution is remaining
- 7 Add 5 ml conc HNO<sub>3</sub> Evaporate until volume is ~1 drop

	RF2	RF7	RF9	RF14	RF19	RF27	RF31	RF32
1	✓	✓	✓	✓	✓	✓	✓	✓
2	✓	✓	✓	✓	✓	✓	✓	✓
3	✓	✓	✓	✓	✓	✓	✓	✓
4								
5								
6	✓	✓	✓	✓	✓	✓	✓	✓
7	✓	✓	✓	✓	✓	✓	✓	✓

**Am Micro-precipitation (SOP 780)**

- 1 Add 1 ml conc HCl to sample Heat for 5 minutes Add 15 ml DI
- 2 Add 0 5 ml lanthanum carrier mix well Add 5 ml HF Mix well
- 3 Allow sample to stand for 15-20 minutes minimum
- 4 Place 25 mm filter membrane in a filter funnel assembly and turn on vacuum Rinse with 1-2 ml alcohol
- 5 Load sample into filter Rinse sample beaker once with 5 ml DI and add to funnel
- 6 After sample has passed through filter Rinse filter with 10-15 ml DI
- 7 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 8 Mount filter on stainless steel planchet with double sided adhesive tape

	RF2	RF7	RF9	RF14	RF19	RF27	RF31	RF32
1	✓	✓	✓	✓	✓	✓	✓	✓
2	✓	✓	✓	✓	✓	✓	✓	✓
3	✓	✓	✓	✓	✓	✓	✓	✓
4	✓	✓	✓	✓	✓	✓	✓	✓
5	✓	✓	✓	✓	✓	✓	✓	✓
6	✓	✓	✓	✓	✓	✓	✓	✓
7	✓	✓	✓	✓	✓	✓	✓	✓
8	✓	✓	✓	✓	✓	✓	✓	✓

Completed 9/28/98



**Pu column (SOP 778)**

- 1 Fill disposable plastic column with 7 cm AG1x8 resin (by resin/DI slurry)
- 2 Place funnel on top of column with Whatman filter paper
- 3 Wet filter paper with 9N HCl
- 4 Condition resin with 50 ml 9 N HCl
- 5 Load sample through filter, rinse 2x with 20 ml 9N HCl (= Am fraction)
- 6 Rinse 2x with 20 ml 9N HCl, discard
- 7 Elute Pu with 20ml 9N HCl + 1.5 ml HI
- 8 Add 1 ml conc HNO<sub>3</sub>, evaporate to dryness

	KF23	KF28	KF29	KF30		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		

**Pu Microprecipitation (SOP 778)**

- 1 Add 1 ml conc HCl, mix well to resolubilize
- 2 Add 14 ml DI, mix well
- 3 Add 1.0 ml lanthanum carrier + 0.5 ml H<sub>2</sub>O<sub>2</sub> mix well
- 4 Add 5 ml 3N HF
- 5 Let sample sit for min 15-20 minutes
- 6 Set up filtration apparatus place 25 mm filter membrane on support screen and lock
- 7 Apply vacuum, rinse filter with 1-2 ml methanol, rinse filter with DI
- 8 Transfer sample, rinsing beaker once with 5 ml DI
- 9 Rinse filter with 10-15 ml DI
- 10 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 11 Mount filter on stainless steel planchet with double sided adhesive tape

	KF23	KF28	KF29	KF30		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		

**Am Methanolic Anion Exchange Column (SOP 780)**

- 1 Mix anion exchange resin with twice the volume of 1N HNO<sub>3</sub>/93% methanol solution (for 160 ml 10ml conc HNO<sub>3</sub>, 150 ml methanol) overnight
- 2 Add 5-10 ml conc HNO<sub>3</sub>,
- 3 Add 100-125 ml DI,
- 4 Add 10 ml Fe carrier
- 5 Transfer to cent tube and ppt with conc NH<sub>4</sub>OH Cent @3500 rpm for 15 min/
- 6 Add conc HNO<sub>3</sub> and evaporate, Repeat
- 7 Add conc HNO<sub>3</sub> until dissolved Add methanol (15 ml of methanol for each 1 ml conc HNO<sub>3</sub>)
- 8 Pour resin slurry into disposable column and let resin settle to 7cm, place a layer of silica or glass beads on top of resin
- 9 Place funnel on top of column with Whatman filter paper, wet filter with 1N HNO<sub>3</sub> /93% methanol solution
- 10 Condition column with 40 ml of 1N HNO<sub>3</sub> /93% methanol solution
- 8 Load the sample onto column through filter, rinse with 25 ml 1N HNO<sub>3</sub> /93% methanol solution
- 9 Remove and discard filter and rinse twice with 25 ml 1N HNO<sub>3</sub> / 93% methanol solution Discard into methanolic waste receptacle
- 10 Strip Am by passing three 20 ml volumes of 8 N HNO<sub>3</sub> through the column allowing each rinse to pass completely before adding next rinse Collect eluate in beaker for Teva column
- 11 Evaporate to dryness

	RF23	RF28	RF29	RF30		
✓	✓	✓	✓			
✓	✓	✓	✓			
9/17/98 →						
✓	✓	✓	✓			
9/21/98 →						
9/23/98 →						
✓	✓	✓	✓			
✓	✓	✓	✓			
9/24/98 →						

**Am Teva Resin (SOP 780)**

- 1 Redissolve sample from above in 10 ml 2M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100ml 15 2 g NH<sub>4</sub>SCN & 0 35 ml 98% formic acid in 100ml DI) Allow sample to sit for 1 hour to ensure dissolution
- 2 Condition a TEVA resin 2 ml column with 5 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid solution
- 3 Transfer sample into TEVA column in two portions using disposable polyethylene transfer pipet Rinse the sample container with 1 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid and transfer to column Repeat rinse and add to column
- 4 Rinse the TEVA column with two 5 ml volumes of 1M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100 ml dissolve 7 6 g NH<sub>4</sub>SCN and 0 35 ml 98% formic acid in 100 ml of DI) Allow first wash to pass completely before adding the second wash (This washes lanthanides from column)
- 5 Strip Am from column with 15 ml of 2N HCl in three 5ml portions allowing each 5 ml to pass completely
- 6 To decompose thiocyanate, add 2 5 ml conc HNO<sub>3</sub> and 7 5 ml conc HCL to the Am solution Swirl gently Evaporate until ~1 drop solution is remaining
- 7 Add 5 ml conc HNO<sub>3</sub> Evaporate until volume is ~1 drop

	RF23	RF28	RF29	RF30		
	9/24/98			→		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	✓	✓	✓	✓		
	9/25/98			→		

**Am Micro-precipitation (SOP 780)**

- 1 Add 1 ml conc HCl to sample Heat for 5 minutes Add 15 ml DI
- 2 Add 0 5 ml lanthanum carrier mix well Add 5 ml HF Mix well
- 3 Allow sample to stand for 15-20 minutes minimum
- 4 Place 25 mm filter membrane in a filter funnel assembly and turn on vacuum Rinse with 1-2 ml alcohol
- 5 Load sample into filter Rinse sample beaker once with 5 ml DI and add to funnel
- 6 After sample has passed through filter Rinse filter with 10-15 ml DI
- 7 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 8 Mount filter on stainless steel planchet with double sided adhesive tape

	RF23	RF28	RF29	RF30		
	9/25/98			→		
	✓	✓	✓	✓		
	-	✓	✓	✓		
	✓		✓	✓		
	✓		✓	✓		
	✓		✓	✓		
	✓		✓	✓		

**Am Methanolic Anion Exchange Column (SOP 780)**

- 1 Mix anion exchange resin with twice the volume of 1N HNO<sub>3</sub>/93% methanol solution (for 160 ml 10ml conc HNO<sub>3</sub>, 150 ml methanol) overnight
- 2 Add 5-10 ml conc HNO<sub>3</sub>,
- 3 Add 100-125 ml DI,
- 4 Add 10 ml Fe carrier
- 5 Transfer to cent tube and ppt with conc NH<sub>4</sub>OH Cent @3500 rpm for 15 min/
- 6 Add conc HNO<sub>3</sub> and evaporate, Repeat
- 7 Add conc HNO<sub>3</sub> until dissolved Add methanol (15 ml of methanol for each 1 ml conc HNO<sub>3</sub>)
- 8 Pour resin slurry into disposable column and let resin settle to 7cm, place a layer of silica or glass beads on top of resin
- 9 Place funnel on top of column with Whatman filter paper, wet filter with 1N HNO<sub>3</sub> /93% methanol solution
- 10 Condition column with 40 ml of 1N HNO<sub>3</sub> /93% methanol solution
- 8 Load the sample onto column through filter, rinse with 25 ml 1N HNO<sub>3</sub> /93% methanol solution
- 9 Remove and discard filter and rinse twice with 25 ml 1N HNO<sub>3</sub> / 93% methanol solution Discard into methanolic waste receptacle
- 10 Strip Am by passing three 20 ml volumes of 8 N HNO<sub>3</sub> through the column allowing each rinse to pass completely before adding next rinse Collect eluate in beaker for Teva column
- 11 Evaporate to dryness

	RF4	RF16	RF1	RF13	RF24
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	9/17/18 →				
	✓	✓	✓	✓	✓
	9/21/18 →				
	9/22/18 →		9/23/18 →		
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
			9/24/18 →		

**Am Teva Resin (SOP 780)**

- 1 Redissolve sample from above in 10 ml 2M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100ml 15 2 g NH<sub>4</sub>SCN & 0 35 ml 98% formic acid in 100ml DI) Allow sample to sit for 1 hour to ensure dissolution
- 2 Condition a TEVA resin 2 ml column with 5 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid solution
- 3 Transfer sample into TEVA column in two portions using disposable polyethylene transfer pipet Rinse the sample container with 1 ml of 2M NH<sub>4</sub>SCN/ 0 1M formic acid and transfer to column Repeat rinse and add to column
- 4 Rinse the TEVA column with two 5 ml volumes of 1M NH<sub>4</sub>SCN/ 0 1M formic acid (for 100 ml dissolve 7 6 g NH<sub>4</sub>SCN and 0 35 ml 98% formic acid in 100 ml of DI) Allow first wash to pass completely before adding the second wash (This washes lanthanides from column)
- 5 Strip Am from column with 15 ml of 2N HCl in three 5ml portions allowing each 5 ml to pass completely
- 6 To decompose thiocyanate, add 2 5 ml conc HNO<sub>3</sub> and 7 5 ml conc HCL to the Am solution Swirl gently Evaporate until ~1 drop solution is remaining
- 7 Add 5 ml conc HNO<sub>3</sub> Evaporate until volume is ~1 drop

	RF4	RF14	RF1	RF13	RF24
	<del>9/22/98</del>		9/24/98		→
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
	<del>9/24/98</del>		9/25/98		→

**Am Micro-precipitation (SOP 780)**

- 1 Add 1 ml conc HCl to sample Heat for 5 minutes Add 15 ml DI
- 2 Add 0 5 ml lanthanum carrier mix well Add 5 ml HF Mix well
- 3 Allow sample to stand for 15-20 minutes minimum
- 4 Place 25 mm filter membrane in a filter funnel assembly and turn on vacuum Rinse with 1-2 ml alcohol
- 5 Load sample into filter Rinse sample beaker once with 5 ml DI and add to funnel
- 6 After sample has passed through filter Rinse filter with 10-15 ml DI
- 7 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 8 Mount filter on stainless steel planchet with double sided adhesive tape

	RF4	RF14	RF1	RF13	RF24
	<del>9/24/98</del>		9/25/98		→
	✓	✓	✓	✓	✓
	✓	✓	✓	✓	✓
			✓	✓	✓
			✓	✓	✓
			✓	✓	✓
			✓	✓	✓

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/02/98 12:09:54 PM

Sample Title : 98826WH  
Spectrum Description :  
Sample Identification : RF1AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/28/98 10:24:34 AM

Live Time : 351938.0 seconds  
Real Time : 351939.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\* \*\*\*\*\*

Detector Name: Q1-DET01

Sample Title: 98826WH

Peak Analysis Performed on: 10/02/98 12:09:53 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	445-	481	469.64	5262.83	1.33E+003	36.48	0.00E+000
2	481-	503	492.12	5455.88	6.00E+001	7.75	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 10:28:14 AM

Sample Title : 98826 TD  
Spectrum Description :  
Sample Identification : RF2Am  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 10/02/98 12:31:20 PM

Live Time : 419025.0 seconds  
Real Time : 419025.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\*  
 P E A K     A N A L Y S I S     R E P O R T     \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q1-DET01  
 Sample Title: 98826 TD  
 Peak Analysis Performed on: 10/12/98 10:28:14 AM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	436-	478	463.94	5213.87	1.20E+003	34.63	0.00E+000
2	481-	503	492.55	5459.58	4.20E+001	6.48	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 9:37:49 AM

Sample Title :  
Spectrum Description :  
Sample Identification : RF3AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/15/98 4:30:50 PM

Live Time : 232713.0 seconds  
Real Time : 232713.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\* \*\*\*\*\*

Detector Name: Q1-DET01

Sample Title:

Peak Analysis Performed on: 10/12/98 9:37:49 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	475	469.82	5264.38	1.04E+002	10.20	0.00E+000
2	477-	503	485.86	5402.13	2.10E+001	4.58	0.00E+000
3	481-	503	492.45	5458.78	1.10E+001	3.32	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*  
\*\*\*\*\*

G A M M A   S P E C T R U M   A N A L Y S I S

\*\*\*\*\*

Report Generated On : 10/28/98 10:37:18 AM

Sample Title : 98826 >20  
Spectrum Description :  
Sample Identification : RF3AM2  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 10/22/98 11:39:33 AM

Live Time : 1009269.0 seconds  
Real Time : 1009271.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q1-DET01

Sample Title: 98826 >20

Peak Analysis Performed on: 10/28/98 10:37:18 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	447-	482	467.37	5243.38	6.96E+002	26.38	0.00E+000
2	483-	505	491.16	5447.69	8.60E+001	9.27	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*>\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*<

Report Generated On : 10/28/98 2:21:51 PM  
Sample Title : 98826 >5  
Spectrum Description :  
Sample Identification : RF4AM  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/23/98 3:10:58 PM  
Live Time : 415469.0 seconds  
Real Time : 415470.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
\*\*\*\*\* P E A K   A N A L Y S I S   R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q3-DET11

Sample Title: 98826 >5

Peak Analysis Performed on: 10/28/98 2:21:51 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	445-	482	470.30	5268.56	1.60E+003	39.97	0.00E+000
2	483-	505	494.83	5480.09	1.29E+003	35.86	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 9:36:52 AM

Sample Title :  
Spectrum Description :  
Sample Identification : RF5AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/15/98 4:30:52 PM

Live Time : 232804.0 seconds  
Real Time : 232805.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \* \*\*\*\*\*

Detector Name: Q1-DET02

Sample Title:

Peak Analysis Performed on: 10/12/98 9:36:52 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	444-	481	468.84	5257.37	7.62E+002	27.60	0.00E+000
2	481-	503	491.73	5451.43	1.10E+001	3.32	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

46

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/28/98 1:23:58 PM

Sample Title : 98826 >0.45 F  
Spectrum Description :  
Sample Identification : RF5AM2  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 10/22/98 11:39:35 AM

Live Time : 1009313.0 seconds  
Real Time : 1009314.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q1-DET02

Sample Title: 98826 >0.45 F

Peak Analysis Performed on: 10/28/98 1:23:58 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	441-	482	468.87	5257.61	3.62E+003	60.14	0.00E+000
2	484-	511	495.28	5481.51	6.90E+001	8.31	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 9:35:01 AM  
Sample Title : 98826 >0.45 C  
Spectrum Description :  
Sample Identification : RF6AM  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/15/98 4:30:54 PM  
Live Time : 232804.0 seconds  
Real Time : 232805.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q1-DET03  
 Sample Title: 98826 >0.45 C  
 Peak Analysis Performed on: 10/12/98 9:35.01 AM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	447-	482	469.77	5265.74	7.02E+002	26.50	0.00E+000
2	481-	503	492.77	5464.73	4.30E+001	6.56	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/28/98 1:37:04 PM  
Sample Title : 98826 >0.45 C  
Spectrum Description :  
Sample Identification : RF6AM2  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/22/98 11:39:37 AM  
Live Time : 1009339.0 seconds  
Real Time : 1009342.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q1-DET03  
 Sample Title: 98826 >0.45 C  
 Peak Analysis Performed on: 10/28/98 1:37:04 PM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	429-	482	469.29	5261.53	3.84E+003	61.98	0.00E+000
2	486-	508	494.52	5479.88	1.91E+002	13.82	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma



\*\*\*\*\*  
 \*\*\*\*\* P E A K    A N A L Y S I S    R E P O R T    \*\*\*\*\*  
 \*    \*\*\*\*\*

Detector Name: Q1-DET02  
 Sample Title: 98826 0.1 U  
 Peak Analysis Performed on: 10/12/98 10:28:40 AM  
     Peak Analysis From Channel: 1  
     Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	466.64	5238.71	1.13E+003	33.66	0.00E+000
2	481-	503	493.11	5463.16	1.80E+001	4.24	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

54

\*\*\*\*\*  
\* G A M M A S P E C T R U M A N A L Y S I S \*  
\*\*\*\*\*

Report Generated On : 10/02/98 12:10:30 PM  
Sample Title : 98826 0.1R  
Spectrum Description :  
Sample Identification : RF8AM  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/28/98 5:02:27 PM  
Live Time : 328068.0 seconds  
Real Time : 328069.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \* \*\*\*\*\*

Detector Name: Q1-DET02  
 Sample Title: 98826 0.1R  
 Peak Analysis Performed on: 10/02/98 12:10:30 PM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	467.89	5249.26	4.55E+002	21.33	0.00E+000
2	481-	503	491.15	5446.56	1.30E+001	3.61	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 10:29:27 AM

Sample Title : 98826 100 U  
Spectrum Description :  
Sample Identification : RF9AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 10/02/98 12:31:23 PM

Live Time : 419138.0 seconds  
Real Time : 419139.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

```
*****  
***** P E A K A N A L Y S I S R E P O R T *****  
* *****
```

Detector Name: Q1-DET03

Sample Title: 98826 100 U

Peak Analysis Performed on: 10/12/98 10:29:27 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	467.39	5245.12	9.61E+002	31.00	0.00E+000
2	481-	503	492.93	5466.10	2.70E+001	5.20	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/02/98 12:10:59 PM

Sample Title : 98826 100R  
Spectrum Description :  
Sample Identification : RF10AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/28/98 10:24:37 AM

Live Time : 351938.0 seconds  
Real Time : 351939.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K    A N A L Y S I S    R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q1-DET03  
 Sample Title: 98826 100R  
 Peak Analysis Performed on: 10/02/98 12:10:59 PM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	469.09	5259.83	1.10E+001	3.32	0.00E+000
2	481-	503	490.00	5440.78	2.00E+000	1.41	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

60

\*\*\*\*\*  
\*\*\*\*\* .            G A M M A   S P E C T R U M     A N A L Y S I S            \*\*\*\*\*  
\*\*\*\*\*

Report Generated On                                : 10/02/98 12:15:49 PM

Sample Title                                        : 98827A WH  
Spectrum Description                              :  
Sample Identification                              : RF13AM  
Sample Type                                        :  
Sample Geometry                                    :

Peak Locate Threshold                            : 3.00  
Peak Locate Range (in channels)                 : 1 - 1024  
Peak Area Range (in channels)                  : 1 - 1024  
Identification Energy Tolerance                 : 0.020 MeV

Sample Size                                        : 1.000E+000 Unit

Sample Taken On                                    :  
Acquisition Started                               : 9/28/98 10:24:39 AM

Live Time                                         : 351938.0 seconds  
Real Time                                         : 351938.0 seconds

Energy Calibration Used Done On                : 8/04/98  
Efficiency Calibration Used Done On            : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\*  
 P E A K     A N A L Y S I S     R E P O R T     \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q1-DET04  
 Sample Title: 98827A WH  
 Peak Analysis Performed on: 10/02/98 12:15:49 PM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	447-	481	469.40	5264.08	1.65E+003	40.56	0.00E+000
2	482-	499	493.07	5465.80	5.90E+001	7.68	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

62

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 8:33:53 AM

Sample Title : 98827A TD  
Spectrum Description :  
Sample Identification : RF14AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 10/07/98 2:26:42 PM

Live Time : 410759.0 seconds  
Real Time : 410760.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q2-DET05

Sample Title: 98827A TD

Peak Analysis Performed on: 10/12/98 8:33:53 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	468.45	5265.83	1.56E+003	39.50	0.00E+000
2	481-	503	494.72	5492.51	3.20E+001	5.66	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

64

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/28/98 9:05:32 AM

Sample Title :  
Spectrum Description :  
Sample Identification : RF15AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/23/98 8:45:51 AM

Live Time : 433155.0 seconds  
Real Time : 433156.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K     A N A L Y S I S     R E P O R T     \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q3-DET09

Sample Title:

Peak Analysis Performed on: 9/28/98 9:05:32 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	483	470.33	5270.08	1.73E+003	41.61	0.00E+000
2	485-	503	493.09	5468.22	1.20E+002	10.95	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

66

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/28/98 2:26:30 PM

Sample Title : 98827A>5  
Spectrum Description :  
Sample Identification : RF16AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 10/23/98 3:10:53 PM

Live Time : 415476.0 seconds  
Real Time : 415476.0 seconds

Energy Calibration Used Done On : 8/06/98  
Efficiency Calibration Used Done On : 8/06/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q3-DET12

Sample Title: 98827A&gt;5

Peak Analysis Performed on: 10/28/98 2:26:30 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	445-	480	468.60	5257.40	1.61E+003	40.06	0.00E+000
2	481-	504	493.53	5470.99	9.24E+002	30.40	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

68

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/28/98 9:06:19 AM

Sample Title :  
Spectrum Description :  
Sample Identification : RF17AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/23/98 8:45:52 AM

Live Time : 433155.0 seconds  
Real Time : 433156.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

v9

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\* \*\*\*\*\*

Detector Name: Q3-DET10

Sample Title:

Peak Analysis Performed on: 9/28/98 9:06:19 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	442-	482	468.73	5266.29	1.87E+003	43.22	0.00E+000
2	484-	503	492.81	5472.63	4.30E+001	6.56	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*  
\*\*\*\*\*

Report Generated On : 9/28/98 9:07:09 AM  
Sample Title :  
Spectrum Description :  
Sample Identification : RF18AM  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/23/98 8:45:55 AM  
Live Time : 433154.0 seconds  
Real Time : 433155.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q3-DET11

Sample Title:

Peak Analysis Performed on: 9/28/98 9:07:09 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	475	468.97	5257.11	1.24E+002	11.14	0.00E+000
2	476-	503	488.73	5427.47	1.43E+002	11.96	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\* \* \* \* \* G A M M A S P E C T R U M A N A L Y S I S \* \* \* \* \*  
\*\*\*\*\*

Report Generated On : 10/28/98 1:37:55 PM  
Sample Title : 98827A >0.45 C  
Spectrum Description :  
Sample Identification : RF18AM2  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/22/98 11:39:39 AM  
Live Time : 1009368.0 seconds  
Real Time : 1009369.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q1-DET04

Sample Title: 98827A >0.45 C

Peak Analysis Performed on: 10/28/98 1:37:55 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	455-	481	470.81	5276.09	2.76E+002	16.61	0.00E+000
2	483-	504	491.26	5450.38	2.40E+002	15.49	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 10:30:00 AM  
Sample Title : 98827A 0.1 U  
Spectrum Description :  
Sample Identification : RF19AM  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/02/98 12:31:25 PM  
Live Time : 419190.0 seconds  
Real Time : 419191.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K    A N A L Y S I S    R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q1-DET04  
 Sample Title: 98827A 0.1 U  
 Peak Analysis Performed on: 10/12/98 10:30:00 AM  
     Peak Analysis From Channel: 1  
     Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	468.06	5252.70	8.21E+002	28.65	0.00E+000
2	481-	503	491.64	5453.65	1.40E+001	3.74	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/28/98 2:13:10 PM

Sample Title : 98827A 0.1 - 0.45  
Spectrum Description :  
Sample Identification : RF20AM2  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 10/23/98 3:10:37 PM

Live Time : 1007579.0 seconds  
Real Time : 1007582.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q2-DET05

Sample Title: 98827A 0.1 - 0.45

Peak Analysis Performed on: 10/28/98 2:13:10 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	453-	474	466.62	5250.02	1.43E+002	11.96	0.00E+000
2	481-	503	493.59	5482.75	1.70E+001	4.12	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 10:36:33 AM

Sample Title : 98827A 100U  
Spectrum Description :  
Sample Identification : RF21AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 10/05/98 10:52:17 AM

Live Time : 414697.0 seconds  
Real Time : 414697.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

79

\*\*\*\*\*  
\*\*\*\*\* P E A K   A N A L Y S I S   R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q3-DET10

Sample Title: 98827A 100U

Peak Analysis Performed on: 10/12/98 10:36:33 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	468.04	5260.39	1.51E+003	38.85	0.00E+000
2	481-	503	493.57	5479.12	2.80E+001	5.29	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

Report Generated On                                : 10/12/98 10:36:57 AM  
  
Sample Title                                        : 98827A 100 R  
Spectrum Description                              :  
Sample Identification                             : RF22AM  
Sample Type                                        :  
Sample Geometry                                  :  
  
Peak Locate Threshold                            : 3.00  
Peak Locate Range (in channels)                : 1 - 1024  
Peak Area Range (in channels)                 : 1 - 1024  
Identification Energy Tolerance                : 0.020 MeV  
  
Sample Size                                        : 1.000E+000 Unit  
  
Sample Taken On                                  :  
Acquisition Started                              : 10/07/98 2:26:53 PM  
  
Live Time                                         : 413052.0 seconds  
Real Time                                         : 413054.0 seconds

Energy Calibration Used Done On                : 8/05/98  
Efficiency Calibration Used Done On            : 8/05/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q3-DET09

Sample Title: 98827A 100 R

Peak Analysis Performed on: 10/12/98 10:36:57 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	469.70	5264.59	6.73E+002	25.94	0.00E+000
2	481-	503	492.70	5464.77	2.30E+001	4.80	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 9:00:03 AM  
Sample Title : 98827A <0.45 filter  
Spectrum Description :  
Sample Identification : RF23AM  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/07/98 2:26:41 PM  
Live Time : 410766.0 seconds  
Real Time : 410766.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

```
*****
*****          P E A K      A N A L Y S I S      R E P O R T          *****
*****
```

Detector Name: Q2-DET08

Sample Title: 98827A <0.45 filter

Peak Analysis Performed on: 10/12/98 9:00:03 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	481	469.27	5264.13	1.26E+003	35.50	0.00E+000
2	481-	503	489.24	5436.78	2.90E+001	5.39	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 8:28:47 AM  
Sample Title : 98827B whole water  
Spectrum Description :  
Sample Identification : RF24AM  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/07/98 2:26:34 PM  
Live Time : 410518.0 seconds  
Real Time : 410518.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q1-DET01

Sample Title: 98827B whole water

Peak Analysis Performed on: 10/12/98 8:28:47 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	467.33	5242.98	7.63E+002	27.62	0.00E+000
2	481-	503	492.58	5459.85	1.90E+001	4.36	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/28/98 9:07:40 AM

Sample Title :  
Spectrum Description :  
Sample Identification : RF26AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/23/98 8:45:56 AM

Live Time : 433155.0 seconds  
Real Time : 433156.0 seconds

Energy Calibration Used Done On : 8/06/98  
Efficiency Calibration Used Done On : 8/06/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q3-DET12

Sample Title.

Peak Analysis Performed on: 9/28/98 9:07:40 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	480	469.10	5261.63	1.31E+003	36.12	0.00E+000
2	481-	503	493.02	5466.60	5.80E+001	7.62	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 8:49:27 AM  
Sample Title : 98827B 0.1 U  
Spectrum Description :  
Sample Identification : RF27AM  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/07/98 2:26:44 PM  
Live Time : 410760.0 seconds  
Real Time : 410760.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K    A N A L Y S I S    R E P O R T    \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q2-DET06  
 Sample Title: 98827B 0.1 U  
 Peak Analysis Performed on: 10/12/98 8:49:27 AM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	425-	478	463.29	5209.10	1.59E+003	39.91	0.00E+000
2	481-	503	491.97	5451.68	3.00E+001	5.48	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

90

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 8:29:30 AM

Sample Title : 98827B 0.1 R  
Spectrum Description :  
Sample Identification : RF28AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 10/07/98 2:26:35 PM

Live Time : 410518.0 seconds  
Real Time : 410518.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q1-DET02  
Sample Title: 98827B 0.1 R  
Peak Analysis Performed on: 10/12/98 8:29:30 AM  
Peak Analysis From Channel: 1  
Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	468.64	5255.68	6.52E+002	25.53	0.00E+000
2	481-	503	491.26	5447.45	3.10E+001	5.57	0.00E+000

M = First peak in a multiplet region  
m = Other peak in a multiplet region  
F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/12/98 8:30:09 AM

Sample Title : 98827B 100U  
Spectrum Description :  
Sample Identification : RF29AM  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 10/07/98 2:26:36 PM

Live Time : 410517.0 seconds  
Real Time : 410518.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \* \*\*\*\*\*

Detector Name: Q1-DET03  
 Sample Title: 98827B 100U  
 Peak Analysis Performed on: 10/12/98 8:30:09 AM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	450-	478	469.55	5263.81	7.86E+002	28.04	0.00E+000
2	481-	503	490.45	5444.65	4.70E+001	6.86	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

94



\*\*\*\*\*  
\* G A M M A S P E C T R U M A N A L Y S I S \*  
\*\*\*\*\*

Report Generated On : 10/28/98 4:33:40 PM  
Sample Title : 98827B 100kDa - 0.45  
Spectrum Description :  
Sample Identification : RF30AM2  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/23/98 3:10:38 PM  
Live Time : 1007615.0 seconds  
Real Time : 1007617.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K    A N A L Y S I S    R E P O R T    \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q2-DET06  
 Sample Title: 98827B 100kDa - 0.45  
 Peak Analysis Performed on: 10/28/98 4:33:40 PM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	447-	482	469.85	5264.57	1.54E+003	39.29	0.00E+000
2	483-	506	492.14	5453.15	5.70E+001	7.55	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/23/98 8:29:58 AM  
Sample Title : 98826WH  
Spectrum Description :  
Sample Identification : RF1PU  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/18/98 4:53:47 PM  
Live Time : 401793.0 seconds  
Real Time : 401794.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q1-DET01  
 Sample Title: 98826WH  
 Peak Analysis Performed on: 9/23/98 8:29:58 AM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	410.42	4754.25	7.70E+001	8.77	0.00E+000
2	435-	465	454.93	5136.49	1.21E+002	11.00	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

100

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/24/98 3:21:06 PM

Sample Title : 98826TD  
Spectrum Description :  
Sample Identification : RF2PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/15/98 4:31:08 PM

Live Time : 409848.0 seconds  
Real Time : 409849.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q3-DET10

Sample Title:

Peak Analysis Performed on: 9/24/98 3:21:06 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	421.98	4865.71	1.54E+003	39.22	0.00E+000
2	435-	465	451.90	5122.04	4.90E+001	7.00	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/15/98 3:57:53 PM  
Sample Title : 98826 >20 F  
Spectrum Description :  
Sample Identification : RF3PU  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/11/98 9:01:35 AM  
Live Time : 370278.0 seconds  
Real Time : 370279.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q1-DET01  
 Sample Title: 98826 >20 F  
 Peak Analysis Performed on: 9/15/98 3:57:53 PM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	413.13	4777.54	1.48E+003	38.52	0.00E+000
2	435-	465	444.02	5042.87	1.29E+002	11.36	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

104

\*\*\*\*\*  
\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*  
\*\*\*\*\*

port Generated On : 10/28/98 2:15:14 PM  
Sample Title : 98826 >20 F  
Spectrum Description :  
Sample Identification : RF3PU2  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/23/98 3:10:45 PM  
Live Time : 1007706.0 seconds  
Real Time : 1007708.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q2-DET07

Sample Title: 98826 >20 F

Peak Analysis Performed on: 10/28/98 2:15:14 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	386-	430	421.63	4856.56	2.69E+003	51.87	0.00E+000
2	435-	460	452.20	5117.25	2.27E+002	15.07	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/23/98 8:30:30 AM

Sample Title : 98826>5 C  
Spectrum Description :  
Sample Identification : RF4PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/18/98 4:53:48 PM

Live Time : 401794.0 seconds  
Real Time : 401795.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q1-DET02

Sample Title: 98826>5 C

Peak Analysis Performed on: 9/23/98 8:30:30 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	421.42	4855.20	1.25E+003	35.37	0.00E+000
2	435-	465	453.20	5124.72	2.25E+003	47.40	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\* \ . \* \* \* \* G A M M A S P E C T R U M A N A L Y S I S \* \* \* \* \*  
\*\*\*\*\*

Report Generated On : 9/15/98 3:58:10 PM

Sample Title : 98826 >0.45 F  
Spectrum Description :  
Sample Identification : RF5PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/11/98 9:01:36 AM

Live Time : 370278.0 seconds  
Real Time : 370279.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q1-DET02

Sample Title: 98826 >0.45 F

Peak Analysis Performed on: 9/15/98 3:58:10 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	421.68	4857.43	1.31E+003	36.19	0.00E+000
2	435-	465	451.36	5109.08	1.40E+001	3.74	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/15/98 3:58:23 PM

Sample Title : 98826 >0.45 C  
Spectrum Description :  
Sample Identification : RF6PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/11/98 9:01:38 AM

Live Time : 370278.0 seconds  
Real Time : 370279.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q1-DET03  
 Sample Title: 98826 >0.45 C  
 Peak Analysis Performed on: 9/15/98 3:58:23 PM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	423.96	4869.26	1.28E+003	35.71	0.00E+000
2	435-	465	454.45	5133.10	1.03E+002	10.15	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

112

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/28/98 2:19:53 PM  
Sample Title : 98827A >5  
Spectrum Description :  
Sample Identification : RF16PU2  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/22/98 11:39:56 AM  
Live Time : 591473.0 seconds  
Real Time : 591474.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K     A N A L Y S I S     R E P O R T     \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q3-DET11  
 Sample Title: 98827A >5  
 Peak Analysis Performed on: 10/28/98 2:19:53 PM  
     Peak Analysis From Channel: 1  
     Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	389-	434	424.41	4872.93	2.34E+003	48.37	0.00E+000
2	434-	465	455.14	5137.87	3.22E+003	56.71	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

114

\*\*\*\*\*  
\* \* \* \* \* G A M M A S P E C T R U M A N A L Y S I S \* \* \* \* \*  
\*\*\*\*\*

Report Generated On : 9/24/98 3:30:37 PM

Sample Title : 98826 01U  
Spectrum Description :  
Sample Identification : RF7PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/15/98 4:31:09 PM

Live Time : 409849.0 seconds  
Real Time : 409850.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K   A N A L Y S I S   R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q3-DET11

Sample Title:

Peak Analysis Performed on: 9/24/98 3:30:37 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
----------	-----------	---------	---------------	--------------	---------------	------------------	------------------

\*\*\*\*\*  
 \*\*\*\*\*'            G A M M A   S P E C T R U M   A N A L Y S I S            \*\*\*\*\*  
 \*\*\*\*\*

Report Generated On                               : 9/23/98 8:31:02 AM

Sample Title                                       : 988260.1R  
 Spectrum Description                             :  
 Sample Identification                           : RF8PU  
 Sample Type                                      :  
 Sample Geometry                                 :

Peak Locate Threshold                          : 3.00  
 Peak Locate Range (in channels)               : 1 - 1024  
 Peak Area Range (in channels)                 : 1 - 1024  
 Identification Energy Tolerance               : 0.020 MeV

Sample Size                                      : 1.000E+000 Unit

Sample Taken On                                 :  
 Acquisition Started                             : 9/18/98 4:53:50 PM

Live Time                                        : 401795.0 seconds  
 Real Time                                        : 401796.0 seconds

Energy Calibration Used Done On               : 8/04/98  
 Efficiency Calibration Used Done On           : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q1-DET03

Sample Title: 988260.1R

Peak Analysis Performed on: 9/23/98 8:31:02 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	424.42	4873.25	1.20E+003	34.64	0.00E+000
2	435-	465	453.47	5124.68	3.80E+001	6.16	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma



\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q3-DET12

Sample Title:

Peak Analysis Performed on: 9/24/98 3:30:57 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	419.62	4837.63	9.10E+001	9.54	0.00E+000
2	435-	465	453.83	5130.84	1.20E+001	3.46	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/23/98 8:31:23 AM

Sample Title : 98826 100R  
Spectrum Description :  
Sample Identification : RF10PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/18/98 4:53:53 PM

Live Time : 401794.0 seconds  
Real Time : 401794.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q1-DET04  
 Sample Title: 98826 100R  
 Peak Analysis Performed on: 9/23/98 8:31:23 AM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	425.14	4886.86	2.71E+003	52.04	0.00E+000
2	435-	465	456.88	5157.37	4.90E+001	7.00	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

122

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q2-DET05

Sample Title: 98827AWH

Peak Analysis Performed on: 9/23/98 8:32:54 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	422.36	4868.18	1.25E+003	35.28	0.00E+000
2	435-	465	454.04	5141.53	7.30E+001	8.54	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

123

\*\*\*\*\*  
\*\*\* " GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/24/98 3:31:12 PM

Sample Title : 98827A TD  
Spectrum Description :  
Sample Identification : RF14PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/15/98 4:30:59 PM

Live Time : 232948.0 seconds  
Real Time : 232948.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q2-DET05

Sample Title:

Peak Analysis Performed on: 9/24/98 3:31:12 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	423.68	4879.58	1.01E+003	31.72	0.00E+000
2	435-	465	454.21	5143.02	1.40E+001	3.74	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

125

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/15/98 3:55:10 PM

Sample Title : 98827A >20 F  
Spectrum Description :  
Sample Identification : RF15PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/13/98 3:00:19 PM

Live Time : 341755.0 seconds  
Real Time : 341755.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q2-DET06  
 Sample Title: 98827A >20 F  
 Peak Analysis Performed on: 9/15/98 3:55:10 PM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	422.49	4863.99	1.08E+003	32.79	0.00E+000
2	435-	465	453.05	5122.52	1.29E+002	11.36	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

127

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/23/98 8:33:17 AM

Sample Title : 98827A >5 C  
Spectrum Description :  
Sample Identification : RF16PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/18/98 4:54:29 PM

Live Time : 401921.0 seconds  
Real Time : 401922.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K    A N A L Y S I S    R E P O R T    \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q2-DET06  
 Sample Title: 98827A >5 C  
 Peak Analysis Performed on: 9/23/98 8:33:17 AM  
                           Peak Analysis From Channel: 1  
                           Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	422.99	4868.23	1.67E+003	40.89	0.00E+000
2	435-	465	454.23	5132.50	2.14E+003	46.24	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/15/98 3:55:37 PM

Sample Title : 98827A >0.45 F  
Spectrum Description :  
Sample Identification : RF17PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/13/98 3:00:21 PM

Live Time : 369805.0 seconds  
Real Time : 369806.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q2-DET07

Sample Title: 98827A >0.45 F

Peak Analysis Performed on: 9/15/98 3:55:36 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	422.72	4865.80	8.26E+002	28.74	0.00E+000
2	435-	465	452.74	5121.87	5.00E+001	7.07	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/15/98 3:56:44 PM  
Sample Title : 98827A >0.45 C  
Spectrum Description :  
Sample Identification : RF18PU  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/13/98 3:00:23 PM  
Live Time : 369806.0 seconds  
Real Time : 369806.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q2-DET08  
 Sample Title: 98827A >0.45 C  
 Peak Analysis Performed on: 9/15/98 3:56:44 PM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	423.27	4866.38	1.32E+003	36.30	0.00E+000
2	435-	465	453.53	5128.00	1.55E+002	12.45	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

133

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/24/98 3:31:34 PM  
Sample Title : 98827 AC 10  
Spectrum Description :  
Sample Identification : RF19PU  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/15/98 4:31:01 PM  
Live Time : 232948.0 seconds  
Real Time : 232948.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q2-DET06

Sample Title:

Peak Analysis Performed on: 9/24/98 3:31:34 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	423.02	4868.45	9.80E+002	31.30	0.00E+000
2	435-	465	453.15	5123.36	1.30E+001	3.61	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

port Generated On : 10/28/98 2:16:12 PM  
Sample Title : 98827A <0.1  
Spectrum Description :  
Sample Identification : RF19PU2  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/23/98 3:10:35 PM  
Live Time : 1007759.0 seconds  
Real Time : 1007763.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

```
*****
***** P E A K   A N A L Y S I S   R E P O R T   *****
*****
```

Detector Name: Q2-DET08

Sample Title: 98827A <0.1

Peak Analysis Performed on: 10/28/98 2:16:11 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	432	423.26	4866.28	3.87E+003	62.20	0.00E+000
2	445-	464	453.54	5128.12	3.50E+001	5.92	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\* G A M M A S P E C T R U M A N A L Y S I S \*  
\*\*\*\*\*

Report Generated On : 9/28/98 9:01:28 AM  
Sample Title : 98827 A O I R  
Spectrum Description :  
Sample Identification : RF20PU  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/23/98 8:45:36 AM  
Live Time : 432965.0 seconds  
Real Time : 432966.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q2-DET05

Sample Title:

Peak Analysis Performed on: 9/28/98 9:01:28 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	424.78	4889.04	2.33E+003	48.22	0.00E+000
2	435-	465	454.92	5149.14	2.60E+001	5.10	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

139

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/28/98 9:02:06 AM

Sample Title : 98827A100J  
Spectrum Description :  
Sample Identification : RF21PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/23/98 8:45:38 AM

Live Time : 432966.0 seconds  
Real Time : 432967.0 seconds

Energy Calibration Used Done On : 8/04/98  
Efficiency Calibration Used Done On : 8/04/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q2-DET06

Sample Title:

Peak Analysis Performed on: 9/28/98 9:02:05 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	423.29	4870.71	1.47E+003	38.33	0.00E+000
2	435-	465	454.43	5134.15	4.20E+001	6.48	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/23/98 8:33:49 AM  
Sample Title : 98827A 100R  
Spectrum Description :  
Sample Identification : RF22PU  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/18/98 4:54:31 PM  
Live Time : 401921.0 seconds  
Real Time : 401921.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q2-DET07  
Sample Title: 98827A 100R  
Peak Analysis Performed on: 9/23/98 8:33:49 AM  
Peak Analysis From Channel: 1  
Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	422.65	4865.27	8.30E+002	28.81	0.00E+000
2	435-	465	454.00	5132.62	1.90E+001	4.36	0.00E+000

M = First peak in a multiplet region  
m = Other peak in a multiplet region  
F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/23/98 8:34:15 AM  
Sample Title : 98827A <0.45 CART FILTERED  
Spectrum Description :  
Sample Identification : RF23PU  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/18/98 4:54:32 PM  
Live Time : 401921.0 seconds  
Real Time : 401921.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q2-DET08

Sample Title: 98827A <0.45 CART FILTERED

Peak Analysis Performed on: 9/23/98 8:34:15 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	424.59	4877.83	1.06E+003	32.51	0.00E+000
2	435-	465	452.38	5118.02	8.00E+000	2.83	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

145

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/23/98 8:34:40 AM

Sample Title : 98827BWH  
Spectrum Description :  
Sample Identification : RF24PU  
Sample Type :  
Sample Geometry :

Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV

Sample Size : 1.000E+000 Unit

Sample Taken On :  
Acquisition Started : 9/18/98 4:54:08 PM

Live Time : 401942.0 seconds  
Real Time : 401944.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q3-DET09  
 Sample Title: 98827BWH  
 Peak Analysis Performed on: 9/23/98 8:34:40 AM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	423.15	4859.34	6.53E+002	25.55	0.00E+000
2	435-	465	456.23	5147.32	2.20E+001	4.69	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

147

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/15/98 3:57:26 PM  
Sample Title : 98827B >20 F  
Spectrum Description :  
Sample Identification : RF26PU  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/13/98 3:00:31 PM  
Live Time : 369891.0 seconds  
Real Time : 369892.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

148

\*\*\*\*\*  
 \*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q3-DET09  
 Sample Title: 98827B >20 F  
 Peak Analysis Performed on: 9/15/98 3:57:26 PM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	424.32	4869.57	1.32E+003	36.39	0.00E+000
2	435-	465	453.92	5127.27	8.00E+001	8.94	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma

149

\*\*\*\*\*  
\*\*\*\*\* G A M M A S P E C T R U M A N A L Y S I S \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 10/28/98 2:16:46 PM  
Sample Title : 98827B <0.1  
Spectrum Description :  
Sample Identification : RF27PU2  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 10/23/98 3:10:55 PM  
Live Time : 1006614.0 seconds  
Real Time : 1006617.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q3-DET09

Sample Title: 98827B <0.1

Peak Analysis Performed on: 10/28/98 2:16:46 PM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	386-	433	424.67	4872.59	3.95E+003	62.87	0.00E+000
2	446-	470	456.36	5148.45	4.20E+001	6.48	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

\*\*\*\*\*  
\*\*\*\*\* GAMMA SPECTRUM ANALYSIS \*\*\*\*\*  
\*\*\*\*\*

Report Generated On : 9/23/98 8:35:02 AM  
Sample Title : 98827B 0.1R  
Spectrum Description :  
Sample Identification : RF28PU  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/18/98 4:54:09 PM  
Live Time : 401943.0 seconds  
Real Time : 401944.0 seconds

Energy Calibration Used Done On : 8/05/98  
Efficiency Calibration Used Done On : 8/05/98

\*\*\*\*\*  
 \*\*\*\*\* P E A K    A N A L Y S I S    R E P O R T    \*\*\*\*\*  
 \*\*\*\*\*

Detector Name: Q3-DET10  
 Sample Title: 98827B 0.1R  
 Peak Analysis Performed on: 9/23/98 8:35:02 AM  
 Peak Analysis From Channel: 1  
 Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	423.23	4876.39	1.70E+003	41.28	0.00E+000
2	435-	465	453.21	5133.32	1.40E+001	3.74	0.00E+000

M = First peak in a multiplet region  
 m = Other peak in a multiplet region  
 F = Fitted singlet

Errors quoted at 1.000 sigma



\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q3-DET11

Sample Title: 98827B 100U

Peak Analysis Performed on: 9/23/98 8:35:19 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	425.50	4882.38	5.66E+002	23.79	0.00E+000
2	435-	465	455.38	5139.92	1.60E+001	4.00	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

155

\*\*\*\*\*  
\*\*\*  
G A M M A S P E C T R U M A N A L Y S I S  
\*\*\*\*\*

Report Generated On : 9/23/98 8:36:23 AM  
Sample Title : 98827B100R  
Spectrum Description :  
Sample Identification : RF30PU  
Sample Type :  
Sample Geometry :  
Peak Locate Threshold : 3.00  
Peak Locate Range (in channels) : 1 - 1024  
Peak Area Range (in channels) : 1 - 1024  
Identification Energy Tolerance : 0.020 MeV  
Sample Size : 1.000E+000 Unit  
Sample Taken On :  
Acquisition Started : 9/18/98 4:54:14 PM  
Live Time : 401941.0 seconds  
Real Time : 401942.0 seconds

Energy Calibration Used Done On : 8/06/98  
Efficiency Calibration Used Done On : 8/06/98

\*\*\*\*\*  
\*\*\*\*\* P E A K A N A L Y S I S R E P O R T \*\*\*\*\*  
\*\*\*\*\*

Detector Name: Q3-DET12

Sample Title: 98827B100R

Peak Analysis Performed on: 9/23/98 8:36:23 AM

Peak Analysis From Channel: 1

Peak Analysis To Channel: 1024

Peak No.	ROI start	ROI end	Peak centroid	Energy (keV)	Net Peak Area	Net Area Uncert.	Continuum Counts
1	383-	433	423.44	4870.38	1.61E+003	40.10	0.00E+000
2	435-	465	453.50	5127.98	2.20E+001	4.69	0.00E+000

M = First peak in a multiplet region

m = Other peak in a multiplet region

F = Fitted singlet

Errors quoted at 1.000 sigma

**Solid Sample Digestion (SOP 773)**

- 1 ~~Combust sample if applicable~~ *Samples in 3N HCl w/ tracers evaporate then to Step 3*
- 2 Add tracers
- 3 In teflon beaker add 20 ml conc HNO<sub>3</sub>, 20 ml conc HCL and 20 ml conc HF
- 4 Ensure solids go in to solution (IF not repeat step 3 ) Evaporate to dryness
- 5 Add 10 ml conc HNO<sub>3</sub> and 5 ml 30% H<sub>2</sub>O<sub>2</sub> Let stand at room temp for a few minutes then evaporate to dryness
- 6 Add 10 ml conc HNO<sub>3</sub>, 100ml DI and 2g H<sub>3</sub>BO<sub>3</sub> Heat to ensure dissolution
- 7 Transfer to 250 ml cent Tube using DI
- 8 Add 1 ml (40mg/ml) Fe carrier
- 9 Ppt Fe(OH)<sub>3</sub> by adding conc NH<sub>4</sub>OH (~50ml)
- 10 Centrifuge at 3500 rpm for 30 minutes
- 11 Dissolve ppt With 3x the vol of ppt with conc HCL
- 12 Add 75 ml 9N HCl
- 13 Add 1 ml saturated NaNO<sub>2</sub> mix thoroughly

	RFH	RF14	RF1	RF13	RF24
9/10/98					→
9/11/98					→
✓	✓	✓	✓	✓	✓
9/12/98					→
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
9/14/98					→

**Pu column (SOP 778)**

- 1 Fill disposable plastic column with 7 cm AG1x8 resin (by resin/DI slurry)
- 2 Place funnel on top of column with Whatman filter paper
- 3 Wet filter paper with 9N HCl
- 4 Condition resin with 50 ml 9 N HCl
- 5 Load sample through filter, rinse 2x with 20 ml 9N HCl (= Am fraction)
- 6 Rinse 2x with 20 ml 9N HCl, discard
- 7 Elute Pu with 20ml 9N HCl + 1.5 ml HI
- 8 Add 1 ml conc HNO<sub>3</sub>, evaporate to dryness

	RF4	CF16	RF1	CF13	RF24
9/16/8	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓

**Pu Microprecipitation (SOP 778)**

- 1 Add 1 ml conc HCl, mix well to resolubilize
- 2 Add 14 ml DI, mix well
- 3 Add 1.0 ml lanthanum carrier + 0.5 ml H<sub>2</sub>O<sub>2</sub>, mix well
- 4 Add 5 ml 3N HF
- 5 Let sample sit for min 15-20 minutes
- 6 Set up filtration apparatus place 25 mm filter membrane on support screen and lock
- 7 Apply vacuum, rinse filter with 1-2 ml methanol, rinse filter with DI
- 8 Transfer sample, rinsing beaker once with 5 ml DI
- 9 Rinse filter with 10-15 ml DI
- 10 Turn off vacuum, use "sharpie" and place dot on outside edge of filter marking which side is up Place filter in pyrex beaker and put in drying oven (90-100°C) for ~ 1min
- 11 Mount filter on stainless steel planchet with double sided adhesive tape

	RF4	CF16	RF1	CF13	RF24
9/16/8	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓

159/159